



Email : sales@laboratoryinstrumentindia.com

Product Name :
Determination Of The Settling Velocity

Product Code :
LBNY-0005-103000243



Description :

Experimental unit to determine the settling velocity of different spheres
Sluice to remove the spheres from the cylinder
2 transparent cylinders
10 spheres of various densities and diameters
Marking of the measuring section
Cover with guide tube to insert the sphere
Stopwatch to measure the sedimentation time
2 areometers to determine the density of the fluids.

Technical Specification :

Technical Data:

2 cylinders
Height: 1330mm each
Sink height: 1000mm each
Inner diameter: 92mm each
Spheres
Aluminium (density: 2,7kg/dm)
2x 10mm diameter
2x 5mm diameter
Polyamide (PA), density: 1,13kg/dm
2x 10mm diameter
Polyoxymethylene (POM), density: 1,41kg/dm

2x 10mm diameter
2x 5mm diameter.



Laboratory Instrument India