



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Determination Of The Settling Velocity

**Product Code :**  
LBNY-0005-103000243



**Description :**

Experimental unit to determine the settling velocity of different spheres  
Sluice to remove the spheres from the cylinder  
2 transparent cylinders  
10 spheres of various densities and diameters  
Marking of the measuring section  
Cover with guide tube to insert the sphere  
Stopwatch to measure the sedimentation time  
2 areometers to determine the density of the fluids.

**Technical Specification :**

**Technical Data:**

2 cylinders  
Height: 1330mm each  
Sink height: 1000mm each  
Inner diameter: 92mm each  
Spheres  
Aluminium (density: 2,7kg/dm)  
2x 10mm diameter  
2x 5mm diameter  
Polyamide (PA), density: 1,13kg/dm  
2x 10mm diameter  
Polyoxymethylene (POM), density: 1,41kg/dm

---

2x 10mm diameter  
2x 5mm diameter.



**Laboratory Instrument India**