



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bernoullis Principle Demonstrator

Product Code :
LBNY-0005-103000248



Description :

Familiarisation with Bernoulli's principle
Axially movable Pitot tube for determining the total pressure at various points within the Venturi nozzle
Venturi nozzle with transparent front panel and measuring points for measuring the static pressures
Water supply using base module or via laboratory supply
6 tube manometers for displaying the static pressures
Flow rate determined by base module
Single tube manometer for displaying the total pressure.

Technical Specification :

Technical Data:

Venturi nozzle
Angle at the outlet: 4
Angle at the inlet: 10,5
A: 84...338mm
Pitot tube
Diameter: 4mm
Movable range: 0...200mm
Measuring ranges
Total pressure: 0...370mm
Static pressure: 0...290mm.



Laboratory Instrument India