



Email : sales@laboratoryinstrumentindia.com

Product Name :
U Steam Generator (3 kW)

Product Code :
LBNY-0002-710007



Description :

U Steam Generator (3 kW)

Technical Specification :

U Steam Generator (3 kW)

It is an electrically heated steam generator (3 kW electric resistance). The resistor has an internal thermostat at 150°C.

This type of generator produces steam for an operating pressure of 0.9 bar, with a max pressure of 2 bar. The operating temperature is 120 °C, which corresponds to saturated steam at 0.9 bar.

It has the necessary accessories for safe use, such as safety valve, pressure switch and level safety.

U Steam Generator (3 kW), U Steam Generator (3 kW) Equipments, U Steam Generator (3 kW) Tools, U Steam Generator (3 kW) Tool Kits, U Steam Generator (3 kW) Manufacturers, U Steam Generator (3 kW) Suppliers from India, China, Kenya



Laboratory Instrument India

