



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydrostatics Trainer

**Product Code :**  
LBNY-0005-103000253



**Description :**

Comprehensive experimental introduction to hydrostatics

1 experimental unit each: measuring the buoyancy force, investigation of the hydrostatic pressure in Liquids, measuring the surface tension, communicating vessels, capillarity

Wide range of accessories included: compressor for generating negative pressures, bottom pressure apparatus, two areometers

Transparent tank for observing the processes.

**Technical Specification :**

**Technical Data:**

Pump

Max. head: 7,6m

Max. flow rate: 9m/h

Power consumption: 250W

Compressor

Pressure at inlet: 240mbar

Pressure at outlet: 2bar

Power: 65W

3 tanks

Supply tank for water: 55L

2 areometers with different measuring ranges

Height 500mm

D=100mm, D=133mm, D=200mm.



**Laboratory Instrument India**