



Email : sales@laboratoryinstrumentindia.com

Product Name :
Marcet Boiler

Product Code :
LBNY-0005-103000256



Description :

Measuring a vapour pressure curve for saturated vapour
Bourdon tube pressure gauge to indicate pressure
Digital temperature display
Boiler with insulating jacket
Temperature limiter and safety valve protect against overpressure in the system.

Technical Specification :

Technical Data:

Boiler, stainless steel: 2L
Safety valve: 20bar
Heater: 2kW
Bourdon tube pressure gauge: -1...24bar
Temperature limiter: 200C
Measuring ranges
Pressure: 0...20bar
Temperature: 0...250C.

Laboratory Instrument India

