



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydrostatic Pressure in Liquids

Product Code :
LBNY-0005-103000257



Description :

Tiltable water tank with fill level scale
Lever arm with different weights
Investigation of the hydrostatic pressure in fluids at rest.

Technical Specification :

Technical Data:

Water tank
Scale: 0...250mm
Content: 0...1,8L
Inclination angle: 0...90
Effective area, max. 75x100mm
Lever arm
Max. length 250mm.



Laboratory Instrument India

