



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Fundamentals of Pressure Measurement

**Product Code :**  
LBNY-0005-103000258



**Description :**

One Bourdon tube pressure gauge each for positive and negative pressure  
Basic experiments for measuring pressure with three different measuring instruments  
Plastic syringe generates test pressures in the millibar range  
U-tube and inclined tube manometer.

**Technical Specification :**

**Technical Data:**

Inclined tube manometer: 0...500mmWC  
Inclined tube manometer: angle 30  
Measuring ranges  
Bourdon tube pressure gauge:  
0...60mbar / -60...0mbar  
U-tube manometer: 0...500mmWC.



**Laboratory Instrument India**

