



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fundamentals of Pressure Measurement

Product Code :
LBNY-0005-103000258



Description :

One Bourdon tube pressure gauge each for positive and negative pressure
Basic experiments for measuring pressure with three different measuring instruments
Plastic syringe generates test pressures in the millibar range
U-tube and inclined tube manometer.

Technical Specification :

Technical Data:
Inclined tube manometer: 0...500mmWC
Inclined tube manometer: angle 30
Measuring ranges
Bourdon tube pressure gauge:
0...60mbar / -60...0mbar
U-tube manometer: 0...500mmWC.



Laboratory Instrument India

