



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Centrifugal Pump Demonstration Unit

**Product Code :**  
LBNY-0005-103000263



**Description :**

A small-scale centrifugal pump demonstration unit, comprising of a water reservoir, the pump, control valves and interconnecting pipe work all mounted on a stainless steel base.

Transparent pump volute for visibility

Supplied with software providing full instructions for setting up, operating, calibrating and performing the teaching exercises. Facilities are provided for logging, processing and displaying data graphically. Full theoretical back-up is included together with a student questions and answers session.

Equipped with electronic measurement sensors for pump head pressure, suction, flow-rate and water temperature.

Capable of being linked to a PC (not supplied) via the interface console (available as an essential accessory).

**Technical Specification :**

Centrifugal Pump Demonstration Unit



**Laboratory Instrument India**