



Email : sales@laboratoryinstrumentindia.com

Product Name :
Series And Parallel Pumps Demonstration Unit

Product Code :
LBNY-0005-103000265



Description :

Experimental unit for series, parallel and individual operation of centrifugal pumps
Measuring of pressure and flow rate
One pump with adjustable speed and one pump with fixed speed
Processing of measured data on a PC software for data acquisition via USB under Windows XP, Windows Vista or Windows 7.

Technical Specification :

Technical Data:
Pumps with motors
Max. head: 10 m
Max. flow rate: 40 L/min with one pump
Pump with adjustable speed: 0...3.300 min⁻¹
Power consumption, both pumps: 650 W
Tank for water: ca. 15L
Pipes and fittings made of PVC.



Laboratory Instrument India

