



Email : sales@laboratoryinstrumentindia.com

Product Name :
Base Service Unit

Product Code :
LBNY-0002-710009



Description :

Base Service Unit

Technical Specification :

Base Service Unit (Common for all available Heat Exchangers)

This unit is common for Heat Exchangers and can work with one or several exchangers.

Heating water in a thermostatic bath.

Pumping the hot water.

Regulation and measurement of the cold water and hot water flows.

Measurement of the inlet and outlet temperatures of the cold water and the hot water.

Measurement of the pressure drop in the exchanger.

Base Service Unit, Base Service Unit Equipments, Base Service Unit Tools, Base Service Unit Tool Kits, Base Service Unit Manufacturers, Base Service Unit Suppliers from India, China, Kenya



Laboratory Instrument India

