



Email : sales@laboratoryinstrumentindia.com

Product Name :
Series Parallel Pumps

Product Code :
LBNY-0005-103000270



Description :

This auxiliary pump is intended to be used in conjunction with the basic Hydraulics Bench. The auxiliary pump is mounted on a support plinth which stands adjacent to the Hydraulics Bench primary pump. This accessory comprises a fixed speed pump assembly and independent discharge manifold interconnected by flexible tubing with quick release connectors.

Technical Specification :

Discharge manifold with flow control valve and pressure gauges
All hoses and fittings supplied for easy connection to hydraulics bench in either series or parallel configuration
Fixed speed pump with similar performance characteristics to pump in Hydraulics Bench
Mounted on floor standing plinth with ON/OFF switch
Educational Software available as an option.



Laboratory Instrument India