




Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydrostatic Pressure

Product Code :
LBNY-0005-103000271



Description :

Distance between balance pan and knife edge: 285mm.
Quadrant, 75mm internal radius, 150mm outer radius, 60mm width.
Balance arm.
Balance pan.
Beam level indicator.
Water tank, Perspex, capacity: 10 liters with level adjustment support screws at the base.
Set of standard weights.
Spirit level.
Adjustable counter-balance.
Computer based learning software.
Overall dimensions: Height: 0.45m, Width: 0.3m, Length: 0.55m.
Services required: Water supply and drainage.

Technical Specification :

Hydrostatic Pressure



Laboratory Instrument India

