




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydrostatic Pressure

**Product Code :**  
LBNY-0005-103000271



**Description :**

Distance between balance pan and knife edge: 285mm.  
Quadrant, 75mm internal radius, 150mm outer radius, 60mm width.  
Balance arm.  
Balance pan.  
Beam level indicator.  
Water tank, Perspex, capacity: 10 liters with level adjustment support screws at the base.  
Set of standard weights.  
Spirit level.  
Adjustable counter-balance.  
Computer based learning software.  
Overall dimensions: Height: 0.45m, Width: 0.3m, Length: 0.55m.  
Services required: Water supply and drainage.

**Technical Specification :**

Hydrostatic Pressure



**Laboratory Instrument India**

