




Email : sales@laboratoryinstrumentindia.com

Product Name :
Flow Meter Demonstration Unit

Product Code :
LBNY-0005-103000272



Description :

Flow Meter Demonstration Unit is a self-contained apparatus to demonstrate the attributes of flow meters utilized in measurement of water flow via pipes or open channels. Metering devices available in different combinations to suit coursework. Reference turbine flowmeter permanently fitted. Service module with reservoir, volumetric measuring tank & pump.

Technical Specification :

Technical Details:

Reservoir capacity: 140 litres
Volumetric tank capacity: 40 litres
Maximum flow rate: 2.3 litres/s
Length of removable test pipe: 750mm
Nominal bore of pipework: 38mm
Pump motor rating: 0.55kW.



Laboratory Instrument India

