



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Precision Pressure Gauge Calibrator

**Product Code :**  
LBNY-0005-103000274



**Description :**

Spindle pump for pressure generation and fine adjustment.  
Piston allows calibration over a wide range of pressures.  
This unit allows pressure gauges to be accurately calibrated up to 300 bar, with an accuracy of +0.015%.  
Priming pump.  
A spirit level (or bubble level) and rotating foot studs are provided to enable the student to level the unit.  
Laboratory calibration certificate.  
Mineral oil is used as the hydraulic fluid.  
Carrying cases for the unit and weights.  
A set of weights, a set of adapters for the pressure gauges, spare seals and a bottle of oil are supplied.

**Technical Specification :**

Precision Pressure Gauge Calibrator



**Laboratory Instrument India**