



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Basic Metacentric Height Apparatus

**Product Code :**  
LBNY-0005-103000277



**Description :**

To determine metacentric height of floating body and height variation with tilt angle. Tilt angle indicated by plumb bob on attached scale. Rectangular pontoon; centre of gravity can be moved side way by moving horizontal jockey weight, and vertically via adjustable vertical weight on mast.

**Technical Specification :**

**Technical Detail:**

Corresponding linear dimension: +90mm  
Max. angle of heel: +13o  
Pontoon dimensions: Length 350mm  
Width: 200mm  
Overall height: 475mm.



**Laboratory Instrument India**