




Email : sales@laboratoryinstrumentindia.com

Product Name :
Laminar Flow Table

Product Code :
LBNY-0005-103000278



Description :

A floor standing unit developed to simulate perfect flow, demonstrate Hele-Shaw principles & accept student made models.

Includes dye injection system and adjust-able knife edge weir for rapid, accurate adjustments of table flow.

A system of pipes, valves and manifolds enables any configuration of sinks and sources to be used.

Working section support and end tanks moulded from glass reinforced plastic (GRP).

Working section made from laminated glass (dimensions of top glass 606mm x 892mm).

Lower glass plate incorporates 8 pre-con-nected sinks/sources in cruciform pattern and viewing graticule.

Includes 3gm blue dye powder to make one litre dye solution and comprehensive in-struction manual with illustration.

Technical Specification :

Laminar Flow Table



Laboratory Instrument India