



Email : sales@laboratoryinstrumentindia.com

Product Name :
Concentric Tube Heat Exchanger

Product Code :
LBNY-0002-710010



Description :

Concentric Tube Heat Exchanger

Technical Specification :

Concentric Tube Heat Exchanger

It has been designed to study the heat transfer between hot water flowing through an internal tube and cold water flowing through the annular area between the internal and the external tubes.

Global energy balance in the heat exchanger and the study of losses.

Exchanger effectiveness determination.

Study of the heat transfer under counter-current and co-current flow conditions.

Flow influence on the heat transfer.

Reynolds number calculation.

Concentric Tube Heat Exchanger, Concentric Tube Heat Exchanger Equipments, Concentric Tube Heat Exchanger Tools, Concentric Tube Heat Exchanger Tool Kits, Concentric Tube Heat Exchanger Manufacturers, Concentric Tube Heat Exchanger Suppliers from India, China, Kenya



Laboratory Instrument India

