



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Osborne Reynolds Apparatus

**Product Code :**  
LBNY-0005-103000279



**Description :**

An easily cleaned dye injection system is incorporated.  
Fluid enters a vertically mounted test section through a carefully profiled bell mouth from a constant head tank.  
Fluid flow is controlled by a needle point globe valve.  
Supply comprises stilling bed, hoses & dye  
Floor standing apparatus to produce classic Osborne Reynolds' experiments.

**Technical Specification :**

Osborne Reynolds Apparatus



**Laboratory Instrument India**