



Email : sales@laboratoryinstrumentindia.com

Product Name :
Fluid Friction Measurements

Product Code :
LBNY-0005-103000282



Description :

A unit for the detailed study of fluid friction head losses which occur when an incompressible fluid flows through pipes, fittings and flow metering devices.

A substantial floor standing tubular steel frame supports test circuits comprising:
4 smooth-bore pipes of different diameters ranging from 4.5mm I.D. to 17.2mm I.D.

artificially roughened pipe

90o bends (large & small radii)

90o elbow

90o mitre

45o elbow, 45 Y, 90 T

Sudden Contraction

Sudden Enlargement

Globe Valve

Perspex Orifice Meter

Perspex Venturi

Gate Valve

Ball Valve

Inline Strainer

38 Tapping Points

Perspex Pipe Section With A Pitot Tube & Static Tapping.

Technical Specification :

Fluid Friction Measurements



Laboratory Instrument India