



Email : sales@laboratoryinstrumentindia.com

Product Name :
Extended Concentric Tube Heat Exchanger

Product Code :
LBNY-0002-710011



Description :

Extended Concentric Tube Heat Exchanger

Technical Specification :

Extended Concentric Tube Heat Exchanger

It has been designed to study the heat transfer between hot water flowing through an internal tube and cold water flowing through the annular area between the internal and the external tubes.

Global energy balance in the heat exchanger and the study of losses.

Exchanger effectiveness determination.

Study of the heat transfer under counter-current and co-current flow conditions.

Flow influence on the heat transfer.

Reynolds number calculation.

Extended Concentric Tube Heat Exchanger, Extended Concentric Tube Heat Exchanger Equipments, Extended Concentric Tube Heat Exchanger Tools, Extended Concentric Tube Heat Exchanger Tool Kits, Extended Concentric Tube Heat Exchanger Manufacturers, Extended Concentric Tube Heat Exchanger Suppliers from India, China, Kenya

Laboratory Instrument India

