



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Fluid Statics And Manometry

**Product Code :**  
LBNY-0005-103000294



**Description :**

Demonstrates the basic principles of hydrostatics and manometry

Single piezometer manometer tube, Scale length 460 mm

Inclined manometer with inclinations of 5°, 30°, 60° and 90° (vertical)

Includes vertical tube with variable cross section, Scale length 460 mm

Includes demonstrations of the following types of manometer:

Enlarged limb-manometer

'U' tube manometer (liquid over liquid), Scale length 460 mm

'U' tube manometer (air over liquid), Scale length 460 mm

Level measurement using Vernier hook and point gauge, Range 0 to 150 mm with 0.1 mm resolution

Allows the effect of friction to be demonstrated when fluid is in motion.

Inverted pressurized U tube manometer, Scale length 460 mm.

**Technical Specification :**

Fluid Statics And Manometry



**Laboratory Instrument India**

