



Email : sales@laboratoryinstrumentindia.com

Product Name :
Compressible Flow Apparatus

Product Code :
LBNY-0005-103000299



Description :

A bench top unit meant to demonstrate & teach the fundamentals of compressible flow to engineering students
Four electronic pressure sensors
Test Sections made from clear acrylic
Complete with convergent-divergent duct capable of achieving Mach-1 velocity at the throat
Advanced torque-vector speed control with electronic torque measurement
Data Logging accessory, complete with educational software and electronic temperature sensor.
Additional test sections available (6 off) complete with bench-top stand
Compressor test accessory available.

Technical Specification :

Features:

Standard unit includes convergent-divergent duct designed to produce Mach 1 velocity at the throat.
Variable speed compressor with accurate electronic speed control
Data Logging option
Electronic pressure sensors.



Laboratory Instrument India

