



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Compressible Flow Apparatus

**Product Code :**  
LBNY-0005-103000299



**Description :**

A bench top unit meant to demonstrate & teach the fundamentals of compressible flow to engineering students  
Four electronic pressure sensors  
Test Sections made from clear acrylic  
Complete with convergent-divergent duct capable of achieving Mach-1 velocity at the throat  
Advanced torque-vector speed control with electronic torque measurement  
Data Logging accessory, complete with educational software and electronic temperature sensor.  
Additional test sections available (6 off) complete with bench-top stand  
Compressor test accessory available.

**Technical Specification :**

**Features:**

Standard unit includes convergent-divergent duct designed to produce Mach 1 velocity at the throat.  
Variable speed compressor with accurate electronic speed control  
Data Logging option  
Electronic pressure sensors.



**Laboratory Instrument India**

