



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Flow Channels

**Product Code :**  
LBNY-0005-103000301



**Description :**

The Flow Channel made up of a clear acrylic working section of large depth to width ratio incorporating undershot and overshot weirs at the inlet and discharge ends respectively.

Models supplied with the channel include broad and sharp crested weirs, large & small diameter cylinders and symmetrical & asymmetrical aerofoils which, in conjunction with the inlet and discharge weirs, permit a varied range of open channel and flow visualisation demonstrations.

A dye injection system incorporated at the inlet to the channel permits flow visualisation in conjunction with a graticule on the rear face of the channel.

Water is fed to the streamlined channel entry via a stilling tank to reduce turbulence. Water discharging from the channel is collected in the volumetric tank of the Hydraulics Bench and returned to the sump for recirculation.

**Technical Specification :**

Six different models for investigation  
Quick release fitting for easy connection to hydraulics bench  
Clear acrylic working section fed from stilling tank  
Dye injection system.



**Laboratory Instrument India**

