



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Propeller Velocity Flowmeter

**Product Code :**  
LBNY-0005-1040006



**Description :**

Suitable for both laboratory and field applications battery powered (rechargeable) and fully portable. Measurements are indicated on the digital display and may be connected to a data logger for later analysis. Propeller Velocity Flowmeter Measurement of velocities in clean conductive fluids within the range 25 to 1500mm/sec (up to 3000mm/sec with high speed probe) Operation in confined spaces with limited intrusive effects.

**Technical Specification :**

Propeller Velocity Flowmeter



**Laboratory Instrument India**