



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rainfall Hydrographs

Product Code :
LBNY-0005-1040007



Description :

A unit designed to obtain catchment rainfall and runoff values as functions of time. A motor driven traversing vessel with seventeen compartments is moved by timer beneath the outlet at a preselected rate to collect the runoff and provide an immediate display of the hydrograph. Comprising a bench- or floor-standing tank with two overhead square pattern spray nozzles supplying water via a flow control valve, flow meter and solenoid valve. The flow range is 0.4 to 4.4 litres/minute. A comprehensive user manual is included in the supply. The tank is 1.2m in length x 0.8m wide x 0.2m deep.

Technical Specification :

Technical Specifications:

Flow meter range: 0.4 - 4.4 l/min
Runoff collector: 17 x 0.5l compartments
Tank dimensions:
Length: 1.2m
Height: 0.2m
Width: 0.6m.

Laboratory Instrument India

