



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tripod Stands

Product Code :
LBNY-0005-1040008



Description :

A mounting tripod constructed from aluminium alloy is supported on three stainless steel rods secured by thumbscrews. The rods are adjustable allowing the stand to be made level.
An adjustable mounting tripod constructed from aluminium alloy, supported by three stainless steel rods secured by thumbscrews. For use with Hook & Point Gauges and Pitot Tubes.
The whole assembly may be adjusted for height on the support rods. A carrier plate mounted on the tripod supports the gauge.
Levelling is assisted by use of a Bullseye spirit level mounted on the top plate.

Technical Specification :

Range: 500mm (nominal)
Diameter of base: 340mm
Overall height: 660mm (without gauge).



Laboratory Instrument India