



Email : sales@laboratoryinstrumentindia.com

Product Name :
Series Liquid Manometers Apparatus

Product Code :
LBNY-0005-10400015



Description :

A range of manometers measuring differential water pressures up to approximately 12.5m H₂O. Scales are graduated in 1mm divisions.

- 1 meter scale differential pressurized water manometer (the air space above the tubes can be pressurized with the pump supplied)
- 1 meter scale differential open water manometer
- 1 meter scale differential water over Mercury manometer
- 500mm scale differential kerosene over water manometer
- 500mm scale differential water over Mercury manometer

Technical Specification :

Series Liquid Manometers Apparatus



Laboratory Instrument India