



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bench For The Study Of The Motor Automation

Product Code :
LBNY-0005-1050005



Description :

The PLC monitors the motors operation and correlates the parameters according to the developed software. The PLCs are the control devices used in automation systems with dedicated software so that they are capable of being programmed to perform corresponding control operations and they are widely used in industrial control and can operate at industry environment conditions.

The performance of the induction asynchronous motor can be tested in two different modes:
Fed by the inverter and PLC
Or fed by the inverter only.

Technical Specification :

The electric machines bench has been designed to satisfy the following basic requirements:
The plan activities
The economic advantage
The number of groups of students who must work simultaneously
The available space of the working area.



Laboratory Instrument India

