



Email : sales@laboratoryinstrumentindia.com

Product Name :
Portable Precision Gauging

Product Code :
LBNY-0005-10500010



Description :

These measuring devices, along with a set of precisely-machined components, are used to complete skills, such as measuring length, taking accurate measurement of inside and outside dimensions, and collecting process data. This Portable Precision Gauging Learning System features calipers, micrometers, and a variety of gauges within a flexible, convenient product for situations where space is too limited for a full-size trainer or when a trainer is needed in multiple environments. Customers can additionally purchase the, which features a snap gauge.

Measurement is a cornerstone of all technical career paths and a major part of quality assurance, which helps to keep product quality high and costs low. Portable Precision Gauging Learning System covers the fundamentals of basic measurement, precision measurement, direct gauging, indirect gauging, and dimensional measurements using both the customary system and the SI metric system all within a mobile, highly durable product.

Features:

World-Class Technical Learning Wherever You Are
High Quality Trainer for On-the-Go Learning
World-Class Multimedia.

Technical Specification :

Portable Precision Gauging



Laboratory Instrument India