



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Coil Vessel Heat Exchanger

**Product Code :**  
LBNY-0002-710016



**Description :**

Coil Vessel Heat Exchanger

**Technical Specification :**

Coil Vessel Heat Exchanger

Allows the study of heat transfer between hot water flowing through a coil and cold water contained in the vessel. The apparatus consists of fabricated shell, inside which coil tubes on outer side are fitted. The present set up is a heat exchanger, with cold water flows through shell.

Valves are provided to control the low rates of hot and cold water. Flow rates of hot and cold water are measured using Rotameters.

A magnetic drive pump is used to circulate the hot water from a re-cycled type water tank, which is fitted with heaters and Digital Temperature Controller.

Coil Vessel Heat Exchanger, Coil Vessel Heat Exchanger Equipments, Coil Vessel Heat Exchanger Tools, Coil Vessel Heat Exchanger Tool Kits, Coil Vessel Heat Exchanger Manufacturers, Coil Vessel Heat Exchanger Suppliers from India, China, Kenya



**Laboratory Instrument India**

