



Email : sales@laboratoryinstrumentindia.com

Product Name :
Slip Ring Three Phase Asynchronous Motor

Product Code :
LBNY-0005-10500036



Description :

Induction motor with both stator and rotor three-phase windings.

Technical Specification :

Technical Features:

Excitation: 92 V, 2 A
Power: 250 W
Voltage: 220/380 V D/Y
Current: 1.65/0.95 A D/Y
Speed: 1350 rpm
Frequency: 50 Hz
Cos : 0.66
Thermal protection.



Laboratory Instrument India