



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Variable Three Phase Transformer

**Product Code :**  
LBNY-0005-10500039



www.laboratoryinstrumentindia.com

**Description :**

The voltage is set by means of a variable auto transformer with rotary knob and the output is floating by means of an isolating transformer with subdivided secondary winding. Variable Three Phase Transformer suitable for fine •step adjustment of three• phase voltages from zero to maximum value. Fitted with mains lamp and motor circuit breaker.

**Technical Specification :**

**Technical Features:**

Rated output: 550 VA  
Secondary phase current: 1.25 A  
Power supply: three phase from mains  
Secondary voltage: 0 to 440 V.



**Laboratory Instrument India**