



Email : sales@laboratoryinstrumentindia.com

Product Name :
Winding Machine For Motors

Product Code :
LBNY-0005-10500042



Description :

Windings:
Coils for motors
Coils:
Maximum diameter of the coils: 180 mm.
Maximum width of the coils: 200 mm.
Wires:
Minimum diameter: 0.1 mm.
Maximum diameter: 2 mm.
Winding modes:
Both automatic and manual winding modes.
Dimensions: 700 x 700 x 1700 (h) mm
Motor:
Three phase motor with brake, 0.75 kW

Technical Specification :

Winding Machine For Motors

Laboratory Instrument India

