



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Winding Machine For Motors And Transformers

**Product Code :**  
LBNY-0005-10500043



**Description :**

Windings:

Coils for motors

Coils for transformers

Coils:

Maximum diameter of the coils: 180 mm.

Maximum width of the coils: 200 mm.

Wires:

Minimum diameter: 0.1 mm.

Maximum diameter: 2 mm.

Winding modes:

Both automatic and manual winding modes.

Wire guide with automatic displacement.

Dimensions: 700 x 700 x 1700 (h) mm

Motor:

Three phase motor with brake, 0.75 kW

**Technical Specification :**

Winding Machine For Motors And Transformers

---

# Laboratory Instrument India

