




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Electrical Machine System

**Product Code :**  
LBNY-0005-10500048



#### **Description :**

##### **Features:**

Modular design provides flexible experimental requirement  
Rotary machine and brake with overheating protection  
Adopting digitized and microprocessor-based measuring instrument to provide high-accuracy measurement  
Brake with constant speed / constant torque function, easy to operate  
Each Module panel height compatible with DIN A4 standard  
Using 4 mm safety sockets and plugs  
Each DC/AC power supply equipped with overload protection  
Training panel uses 5 mm isolation bake lite, printed component symbol, value and function, easy to connect  
Fully protected system safe to connect various kind of machines  
Providing powerful computer measuring software for saving graphic file, drawing and printing characteristic curves  
Connecting to PC, measuring and drawing characteristic curve available  
300W-grade designed equipment suitable for learning the theory and characteristics of electrical machines  
Stand-alone machine design equipped with two shaft ends and aluminum alloy base for coupling to other machines  
For the sake of safety, the system normally operates at three-phase 220V. Different line voltages can be adjusted by system transformer.

#### **Technical Specification :**

Electrical Machine System



**Laboratory Instrument India**