



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Introduction To Refrigeration Unit

**Product Code :**  
LBNY-0002-710020



**Description :**

Introduction To Refrigeration Unit

**Technical Specification :**

Introduction To Refrigeration Unit

With simple and clear experiments, the Introduction to Refrigeration Unit, allows the familiarization with the operation and handling of a refrigeration system and offers a basic introduction to refrigeration.

Trainer is designed with simple and functional concepts to visualize the main parts of compression refrigeration cycle. The total set-up is housed in well-designed rigid structure.

Experiments:

To calculate refrigerant capacity.

To study the vapor compression refrigerator cycle.

To calculate COP.

Utilities Required:

Stabilized power supply 15 Amp 220 V AC

Water Supply & Drain.

Introduction To Refrigeration Unit, Introduction To Refrigeration Unit Equipments, Introduction To Refrigeration Unit Tools, Introduction To Refrigeration Unit Tool Kits, Introduction To Refrigeration Unit Manufacturers, Introduction To Refrigeration Unit Suppliers from India, China, Kenya



**Laboratory Instrument India**