



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Digital-Analog Training System

**Product Code :**  
LBNY-0005-10600018



**Description :**

This training system is used as essential part of circuit trouble shooting course, analog circuit experimental work and basic logical circuit training procedure. Design of this Digital-Analog Training System has been developed to meet various research and development laboratories where digital and analog training programs are held. We are a reliable supplier and manufacturer of Digital-Analog Training System. It consists of instruments and sections like potentiometer, BAR display unit, LED display screen and adapter display module to name a few. Besides these, this system comprises of in built function generator, easy to dismantle breadboard and DC power supply system.

**Features:**

Application of latest technology  
Easy to handle  
Precise design.

**Technical Specification :**

Digital-Analog Training System



**Laboratory Instrument India**

