



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Logic Trainer

**Product Code :**  
LBNY-0005-10600020



**Description :**

It is equipped with various kinds of basic logic gates, debounced logical switches, LED indicators, DC power supply with short circuit protection, pulse generator and solderless bread board. The design is easy to operate and easy to understand. Logic Trainer: The logic trainer is designed for the logic beginners to enhance the comprehension of basic logical theory.

**Features:**

- Pulse generator
- Eight LED indicators
- Debounced logic switches
- Basic logic gates
- DC power supply +5V, -5V and  $\pm 15V$
- Solderless breadboard
- Console style case.

**Technical Specification :**

Logic Trainer

---

# Laboratory Instrument India

