



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cheese Vat And Cheese Making Accessories

**Product Code :**  
LBNY-0005-1070004



**Description :**

Manufactured wholly in stainless steel and bench mounted, the vat is hygienic in use. Water is circulated around the vat jacket at pre-selectable temperatures by an internal pump. It has a removable paddle agitator of stainless steel, driven by a variable speed motor.

Digital display of temperature of vat contents and water temperature

Agitator stop/start and variable speed control

The control system incorporates

Digital display of pH of vat contents

Integrated data logger with USB output

A specially designed perforated strainer fits into the vat outlet to facilitate efficient draining.

Automatic control of agitator speed, temperatures and times for up to six steps

The pH and temperature probes are suspended in the vat contents from a bracket on the rim of the vat.

**Technical Specification :**

Cheese Vat And Cheese Making Accessories



**Laboratory Instrument India**

