



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stirred Refrigerated Tank

Product Code :
LBNY-0005-10700014



Description :

Inner tank and outer casing made of stainless steel
Inner trap for emptying
Stainless steel measuring stick for the measuring of the milk level
DN50 threaded outlet connector with plug
Total primary surface evaporator with maximum yield
Non polluting polyurethane foam

Cooling System:

Compact design with space saving installation
Non icing cooling even with low filling level
R404A refrigerating fluid for all models
Direct expansion.

Technical Specification :

Electric System:

Computerised milk controller system for the management of all operations:
Digital temperature display
Stirring program
Milk temperature check
Electrical feeding cables
230V 50HZ single phase feeding.



Laboratory Instrument India