



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Maintenance Kit containing Leakage Detector

**Product Code :**  
LBNY-0002-710022



**Description :**

Maintenance Kit containing Leakage Detector

**Technical Specification :**

Maintenance Kit containing Leakage Detector.

The kit additionally has a LED Flashlight and UV protected glasses for operator safety

The Maintenance Kit containing Leakage Detector, includes a portable detector sensitive to concentrations of all commercial refrigerants with the objective of locating leaks in refrigeration systems. Refrigerants detectable by the kit include HFC, HCFC and CFC.

Detects all HFC (134a) CFC refrigerants; Protect the environment!

Audio 'siren' and visual LED light alert when leak is found; Sensitivity adjustment dial; Flexible probe allows you to reach into awkward corners and tight spaces

Highly precise for accurate readings time after time; Re-set function enables you to set the meter to ignore ambient concentrations of halogen

Reads in just 6 seconds; Battery life indicator

Check your entire refrigerating system including the fittings, house to line couplings, controllers, service parts, the areas around the points of attachment, hold downs on lines and components.

Maintenance Kit containing Leakage Detector, Maintenance Kit containing Leakage Detector Equipments, Maintenance Kit containing Leakage Detector Tools, Maintenance Kit containing Leakage Detector Tool Kits, Maintenance Kit containing Leakage Detector Manufacturers, Maintenance Kit containing Leakage Detector Suppliers from India, China, Kenya



**Laboratory Instrument India**