



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Butter Production Plant

**Product Code :**  
LBNY-0005-10700021



**Description :**

This pilot plant is designed to produce butter by using the cream output by pilot plant. It must be fed with the cream separated from at least 200 litres of milk.

**Technical Specification :**

**Technical Specifications:**

Tank of stainless steel, with capacity of 50 litres

Rotor of stainless steel, equipped with motor

Proximity contact for lid

Framework of stainless steel

Switchboard of stainless steel, with protection degree IP55 including a differential circuit breaker and a schematic diagram of the plant

Max. feeding : 19 l of cream per working cycle (higher flow rates on demand)

Emergency button

Power supply: 230 Vac 50 Hz single-phase - 1 kVA

Dimensions: 500×600×1900 mm

Weight: 150 kg.

---

# Laboratory Instrument India

