



Email : sales@laboratoryinstrumentindia.com

Product Name :
Yoghurt Production Unit

Product Code :
LBNY-0005-1070022



Description :

This plant is designed to produce yoghurt from whole milk or from semi-skimmed milk with addition of enzymes.

Technical Specification :

Technical Specifications:

Wheeled framework of stainless steel

Centrifugal feed pump of stainless steel being also used for recycling

Jacketed fermenter of stainless steel, with capacity of 150 .

Max. feeding: 100 l of milk per working cycle (higher flow rates on demand)

Agitator of stainless steel

Motor variable-speed drive for agitator

Temperature control valve

Electronic temperature control

Switchboard of stainless steel with protection degree complying with EC standards and including a differential circuit breaker and a schematic diagram of the plant

Emergency button

Dimensions: 1500 × 900 × 1900 mm

Weight: 230 kg.

Laboratory Instrument India

