



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tall Form Spray Dryer And Chiller

Product Code :
LBNY-0005-10700023



Description :

Compact, mobile design
Easily dismantled for inspection and cleaning purposes
Tall Form chamber straight side to diameter ratio of 3:1, conical discharge section
Advanced 2-fluid nozzle atomisation system
Data logging facility for all key operational parameters
Control & measurement of relative humidity
Spark arrester on input and enhanced safety extraction fan
Powder discharges from two points, chamber and cyclone
Low noise levels
Accessory sets are available to quickly change from drying to spray chilling
Individual inlet and exhaust fans with variable speed control.
Supplied with optional nozzle configurations for counter- and co-current atomising.

Technical Specification :

Tall Form Spray Dryer And Chiller



Laboratory Instrument India

