



Email : sales@laboratoryinstrumentindia.com

Product Name :
Spray Dryer

Product Code :
LBNY-0005-10700024



Description :

The unit is supplied complete; ready for operation.

Materials of construction:

Drier components: glass with gasket free ground glass flanges

Cabinet: all steel, coated with hard chemical resistant finish

Control panel includes:

Rotary knob to control pump speed

Digital controller and display of inlet temperature

On/off switches for the blower, heater and mains supply

Digital display of outlet temperature

Pilot lights positioned in the panel indicate specific operations to be illustrated.

Technical Specification :

Spray Dryer



Laboratory Instrument India

