



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cheese Vat And Cheese Making

**Product Code :**  
LBNY-0005-10700028



**Description :**

Water heater fitted in console  
Temperature and pH can be monitored by microcomputer using an appropriate interface  
Ten litre capacity stainless steel jacketed vat  
Bench mounting with small footprint  
Simple controls  
Range of cheese making accessories  
Separate console controls vat temperature and paddle speed.

**Technical Specification :**

Cheese Vat And Cheese Making



**Laboratory Instrument India**